

2026

PRODUCT CATALOGUE

WWW.PEKYBIO.COM



“ Let’ s Do Together To Make Lab Research Easier! ”

PEKYBIO® is a globally recognized leader in designing and manufacturing premium laboratory consumables.

Our core product lines include:

- Sample Handling - Bottle Filtration System, Microcentrifuge Tube, Centrifuge Tube, Transport Tube, Micro-Sample Transfer Scoop, Micro Plastic Tube, Cryogenic Storage Tube, PCR Tube & PCR Plate, Deep Well Plate, Silicone Sealing Mat, Aluminum Foil Sealing Film.
- Liquid Handling - Serological Pipette, Aspirating Pipette, Conductive Automation Tip, Universal Pipette Tip, LCTS Pipette Tip, Pasteur Pipette, Exact Volume Transfer Pipette.
- Cell Culture - Cell Culture Plate, Cell Culture Flask, Cell Culture Dish, Petri Dish, Erlenmeyer Flask.
- Sample Storage - PET/PETG Media Bottle, HDPE/PP Reagent Bottle.

All products are manufactured in our ISO-certified facilities with rigorous quality control from raw materials to finished goods.

PEKYBIO® offers factory-trained customer service and technical support.

All PEKYBIO® products are available through leading laboratory supply distributors worldwide.

Thank you for your interest in PEKYBIO® laboratory consumables. We’re pleased to introduce our life science products engineered to support your research. Wherever you are and whatever you need, we’re honored to be part of your success.

For more information, call PEKYBIO at +86 571 87566925 or visit www.pekybio.com.





ZHEJIANG PEKY BIOTECH CO., LTD.



LIFE SCIENCES[®]
PEKYBIO

CONTENTS

01

SAMPLE HANDLING

Bottle Filtration System	7-8
Microcentrifuge Tube	9-10
Centrifuge Tube	11-13
Transport Tube	14-15
Micro-Sample Transfer Scoop	15
Micro Plastic Tube	16-17
Cryogenic Storage Tube	18-19
PCR Tube & PCR Plate	20-22
Deep Well Plate	23-24
Silicone Sealing Mat	24
Aluminum Foil Sealing Film	24

02

LIQUID HANDLING

Serological Pipette	26-27
Aspirating Pipette	28
Conductive Automation Tip	29
Universal Pipette Tip	30-36
LCTS Pipette Tip	37-40
Pasteur Pipette	41-42
Exact Volume Transfer Pipette	42

03

CELL CULTURE

Cell Culture Plate	44-45
Cell Culture Flask	46-47
Cell Culture Dish	48-49
Petri Dish	49-50
Erlenmeyer Flask	51-52

04


SAMPLE STORAGE

Media Bottle	54-55
Reagent Bottle	56-60



01 SAMPLE HANDLING

Bottle Filtration System/ 7-8
Microcentrifuge Tube/ 9-10
Centrifuge Tube/ 11-13
Transport Tube/ 14-15
Micro-Sample Transfer Scoop/ 15
Micro Plastic Tube/ 16-17
Cryogenic Storage Tube/ 18-19
PCR Tube & PCR Plate/ 20-22
Deep Well Plate/ 23-24
Silicone Sealing Mat/ 24
Aluminum Foil Sealing Film/ 24



BOTTLE FILTRATION SYSTEM

PEKYBIO Bottle Filtration Systems are engineered for efficient particle separation in laboratory liquid samples. These systems feature USP Class VI polystyrene construction for superior clarity and chemical resistance. The low-protein-binding, low-extractables membrane delivers optimal performance for sterile filtration of cell culture media, serum, additives and buffers. With high flow rates and exceptional flux, the systems significantly reduce filtration time while maximizing efficiency.



MULTIPLE MEMBRANE SOLUTIONS

MOLDED-IN GRADUATIONS



GL45 STANDARD

INDIVIDUAL PACKAGING



- ◆ Manufactured from premium polystyrene.
- ◆ PES membrane features low protein binding.
- ◆ Fast flow rates and high throughput enhance filtration efficiency.
- ◆ Clear graduations are molded on both filter cup and receiver bottle.
- ◆ Secure component integration ensures firmness and leak-proof sealing.
- ◆ Sterile individual packaging prevents contamination during opening.

The most efficient single-use, disposable process for sterile media preparation.



01 SAMPLE HANDLING

Bottle Filtration System- 150, 250, 500, 1000 mL



Cat. No.	Volume	Membrane	Pore Size	Sterile	Packaging	
PK00775	150 mL	PES	0.22 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK00776	150 mL	PES	0.45 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK00560	250 mL	PES	0.22 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK00561	250 mL	PES	0.45 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK00927	500 mL	PES	0.22 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK00928	500 mL	PES	0.45 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK01080	1000 mL	PES	0.22 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK01081	1000 mL	PES	0.45 µm	YES	Individual Packaging	1piece/bag. 12bags/carton

For other membranes, please contact us if needed.

Bottle Top Vacuum Filter

Receiver Bottle with Closed Cap



Bottle Top Vacuum Filter- 150, 250, 500, 1000 mL



Cat. No.	Volume	Membrane	Pore Size	Sterile	Packaging	
PK00951	150 mL	PES	0.22 µm	YES	Individual Packaging	1piece/bag. 24bags/carton
PK00952	150 mL	PES	0.45 µm	YES	Individual Packaging	1piece/bag. 24bags/carton
PK00953	250 mL	PES	0.22 µm	YES	Individual Packaging	1piece/bag. 24bags/carton
PK00954	250 mL	PES	0.45 µm	YES	Individual Packaging	1piece/bag. 24bags/carton
PK00959	500 mL	PES	0.22 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK00960	500 mL	PES	0.45 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK00961	1000 mL	PES	0.22 µm	YES	Individual Packaging	1piece/bag. 12bags/carton
PK00962	1000 mL	PES	0.45 µm	YES	Individual Packaging	1piece/bag. 12bags/carton

Receiver Bottle with Closed Cap- 150, 250, 500, 1000 mL



Cat. No.	Volume	Cap Size	Cap Color	Sterile	Packaging	
PK01136	150 mL	45 mm	Orange	YES	Individual Packaging	1piece/bag. 24bags/carton
PK01137	250 mL	45 mm	Orange	YES	Individual Packaging	1piece/bag. 24bags/carton
PK01138	500 mL	45 mm	Orange	YES	Individual Packaging	1piece/bag. 12bags/carton
PK01139	1000 mL	45 mm	Orange	YES	Individual Packaging	1piece/bag. 12bags/carton




MICROCENTRIFUGE TUBE

PEKYBIO Microcentrifuge Tubes are crafted from transparent polypropylene, engineered for compatibility with centrifuges and widely employed in molecular biology, clinical chemistry, and biochemical research. Designed to meet routine sample handling requirements, these tubes feature a secure seal that prevents leakage while ensuring effortless opening. Available in four standard volumes: 0.5mL, 1.5mL, 2.0mL, and 5.0mL.



OPTIMAL WRITING PERFORMANCE

MOLDED-IN GRADUATIONS



OPTIMAL SEALING PERFORMANCE

RCF TO 20000xg



- Engineered for single-handed operation with thumb-operated cap.
- Optically clear walls feature a frosted labeling surface for precise sample identification.
- Certified RNase/DNase-free and Non-pyrogen.
- Withstand up to 20,000xg relative centrifugal force.



0.5mL



1.5mL



2.0mL



5.0mL

01 SAMPLE HANDLING


Microcentrifuge Tube- 0.5, 1.5, 2.0, 5.0 mL



Cat. No.	Volume	Color	Type of Cap	Type of Bottom	Sterile	Packaging
PK00123	0.5 mL	Clear	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00124	1.5 mL	Clear	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00124G	1.5 mL	Green	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00124P	1.5 mL	Purple	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00124B	1.5 mL	Blue	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00124R	1.5 mL	Red	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00124Y	1.5 mL	Yellow	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00124O	1.5 mL	Orange	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00124A	1.5 mL	Amber	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00125	2.0 mL	Clear	Attached Flat Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00126	5.0 mL	Clear	Attached Flat Cap	Conical Bottom	NO	100pieces/bag. 20bags/carton

Microcentrifuge Tube with Safe-Lock System

- Special Safe-Lock to ensure complete leak-proof.
- Frosted surface for sample identification.



Microcentrifuge Tube with Safe-Lock System- 1.5, 2.0 ,5.0mL



Cat. No.	Volume	Color	Type of Cap	Type of Bottom	Sterile	Packaging
PK00528	1.5 mL	Clear	Attached Lock Cap	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK00529	2.0 mL	Clear	Attached Lock Cap	Round Bottom	NO	500pieces/bag. 10bags/carton
PK01169	5.0 mL	Clear	Attached Lock Cap	Conical Bottom	NO	100pieces/bag. 20bags/carton



CENTRIFUGE TUBE

PEKYBIO Centrifuge Tubes meet routine sample handling demands across molecular biology, clinical chemistry, and environmental labs. Centrifugation effectively separates particulates from liquids, enabling efficient supernatant removal. Engineered from clear polypropylene for superior chemical resistance, our tubes feature amber options for light-sensitive samples and a ICP/MS variant to eliminate trace contamination. Available in standard volumes: 15mL and 50mL.



PREMIUM RAW MATERIAL

**ANTI-SLIP
ERGONOMIC CAP**



**OPTIMAL WRITING
PERFORMANCE**

CLEAR GRADUATIONS




- ◆ Certified RNase/DNase-free and Non-pyrogen.
- ◆ Anti-slip ergonomic cap ensures leak-proof sealing.
- ◆ Clear graduations for precise reading.
- ◆ Meet bioanalytical-grade standards.
- ◆ Large printed white surface for sample identification and recording.

15 mL



50 mL



**50 mL
Self-Standing**



01 SAMPLE HANDLING

Centrifuge Tube-15, 50 mL



Cat. No.	Volume	Color of Cap	Type of Cap	Color of Tube	Type of Bottom	Sterile	Packaging
Bag/25pcs Packaging							
PK00127	15 mL	Orange	Screw Cap	Natural	Conical Bottom	YES	25pieces/bag. 20bags/carton
PK00128	50 mL	Orange	Screw Cap	Natural	Conical Bottom	YES	25pieces/bag. 20bags/carton
PK00129	50 mL	Orange	Screw Cap	Natural	Self Standing	YES	25pieces/bag. 20bags/carton
PK00370	15 mL	Orange	Screw Cap	Natural	Conical Bottom	NO	25pieces/bag. 20bags/carton
PK00371	50 mL	Orange	Screw Cap	Natural	Conical Bottom	NO	25pieces/bag. 20bags/carton
PK00372	50 mL	Orange	Screw Cap	Natural	Self Standing	NO	25pieces/bag. 20bags/carton
PK01257	15 mL	Orange	Screw Cap	Amber	Conical Bottom	YES	25pieces/bag. 20bags/carton
PK01258	50 mL	Orange	Screw Cap	Amber	Conical Bottom	YES	25pieces/bag. 20bags/carton
PK01259	50 mL	Orange	Screw Cap	Amber	Self Standing	YES	25pieces/bag. 20bags/carton
PK01260	15 mL	Orange	Screw Cap	Amber	Conical Bottom	NO	25pieces/bag. 20bags/carton
PK01261	50 mL	Orange	Screw Cap	Amber	Conical Bottom	NO	25pieces/bag. 20bags/carton
PK01262	50 mL	Orange	Screw Cap	Amber	Self Standing	NO	25pieces/bag. 20bags/carton
Bag/500pcs Packaging							
PK01124	15 mL	Orange	Screw Cap	Natural	Conical Bottom	NO	500pieces/bag. 1bag/carton
PK01125	50 mL	Orange	Screw Cap	Natural	Conical Bottom	NO	500pieces/bag. 1bag/carton
PK01126	50 mL	Orange	Screw Cap	Natural	Self Standing	NO	500pieces/bag. 1bag/carton
Paperboard Rack Packaging							
PK01120	15 mL	Orange	Screw Cap	Natural	Conical Bottom	YES	25pieces/rack. 20racks/carton
PK01121	50 mL	Orange	Screw Cap	Natural	Conical Bottom	YES	25pieces/rack. 20racks/carton
PK01263	15 mL	Orange	Screw Cap	Amber	Conical Bottom	YES	25pieces/rack. 20racks/carton
PK01264	50 mL	Orange	Screw Cap	Amber	Conical Bottom	YES	25pieces/rack. 20racks/carton

ICP/MS Centrifuge Tube

PEKYBIO ICP/MS Centrifuge Tubes are manufactured from high-purity polypropylene via a validated manufacturing process to minimize elemental background contamination during production and handling. Enhanced material purity and process control ensure consistent performance in trace metal analysis applications including ICP-MS, ICP-OES, environmental analysis, battery material characterization and pharmaceutical impurity testing.

- ◆ High-purity polypropylene.
- ◆ Trace elements below 1 ppb (10^{-9}).
- ◆ Tested for over 30 common trace elements.
- ◆ Certified RNase/DNase-free and Non-pyrogen.
- ◆ Sterilized by irradiation, $SAL=10^{-6}$
- ◆ Large printed white surface for sample identification and recording.



15 mL



50 mL



50 mL Self-Standing




ICP/MS Centrifuge Tube- 15, 50 mL



Cat. No.	Volume	Color of Cap	Type of Cap	Type of Bottom	Sterile	Packaging
PK01209	15 mL	Natural	Screw Cap	Conical Bottom	YES	25pieces/bag. 20bags/carton
PK01210	50 mL	Natural	Screw Cap	Conical Bottom	YES	25pieces/bag. 20bags/carton
PK01211	50 mL	Natural	Screw Cap	Self Standing	YES	25pieces/bag. 20bags/carton
PK01212	15 mL	Natural	Screw Cap	Conical Bottom	NO	25pieces/bag. 20bags/carton
PK01213	50 mL	Natural	Screw Cap	Conical Bottom	NO	25pieces/bag. 20bags/carton
PK01214	50 mL	Natural	Screw Cap	Self Standing	NO	25pieces/bag. 20bags/carton

PEKYBIO Centrifuge Tube Production Line



01 SAMPLE HANDLING



TRANSPORT TUBE

PEKYBIO Transport Tubes are engineered for seamless integration from sample collection to laboratory analysis. Widely adopted in water quality monitoring and environmental monitoring workflows, they ensure sample integrity from field deployment to final analytical results.

- ◆ Self-standing design with conical bottom for efficient sample recovery.
- ◆ Meet air transport compliance standards.
- ◆ Optically clear tube with molded-in graduations enable clear visualization of reaction component sedimentation.
- ◆ Superior sealing performance prevents both liquid leakage and sample contamination.




Natural Cap



Blue Cap



5 mL



10 mL

Transport Tube- 5, 10 mL



Cat. No.	Volume	Color of Cap	Type of Cap	Color of Tube	Sterile	Packaging
Cap and Tube Pre-assembly						
PK00144	5 mL	Natural	Screw Cap	Natural	NO	100pieces/bag. 10bags/carton
PK00154	5 mL	Natural	Screw Cap	Natural	YES	100pieces/bag. 10bags/carton
PK00145	5 mL	Blue	Screw Cap	Natural	NO	100pieces/bag. 10bags/carton
PK00155	5 mL	Blue	Screw Cap	Natural	YES	100pieces/bag. 10bags/carton
PK00146	10 mL	Natural	Screw Cap	Natural	NO	100pieces/bag. 10bags/carton
PK00156	10 mL	Natural	Screw Cap	Natural	YES	100pieces/bag. 10bags/carton
PK00147	10 mL	Blue	Screw Cap	Natural	NO	100pieces/bag. 10bags/carton
PK00157	10 mL	Blue	Screw Cap	Natural	YES	100pieces/bag. 10bags/carton
Cap and Tube Separation						
PK00408	5 mL	Natural	Screw Cap	Natural	NO	200pieces/bag. 5bags/carton
PK00412	5 mL	Natural	Screw Cap	Natural	YES	200pieces/bag. 5bags/carton
PK00409	5 mL	Blue	Screw Cap	Natural	NO	200pieces/bag. 5bags/carton
PK00413	5 mL	Blue	Screw Cap	Natural	YES	200pieces/bag. 5bags/carton
PK00410	10 mL	Natural	Screw Cap	Natural	NO	100pieces/bag. 10bags/carton
PK00414	10 mL	Natural	Screw Cap	Natural	YES	100pieces/bag. 10bags/carton
PK00411	10 mL	Blue	Screw Cap	Natural	NO	100pieces/bag. 10bags/carton
PK00415	10 mL	Blue	Screw Cap	Natural	YES	100pieces/bag. 10bags/carton

MICRO-SAMPLE TRANSFER SCOOP

PEKYBIO Micro-Sample Transfer Scoops are designed specifically for the transfer of powders in small quantities, making them suitable for a wide range of laboratory applications.

- Antistatic polypropylene minimizes sample loss during powder transfer.
- Micro-scoops deliver 2-12 mg precision for consistent small-volume powder sample handling.
- High-purity polypropylene ensures chemical resistance and autoclavability.
- Disposable design eliminates cross-contamination risk.



Micro-Sample Transfer Scoop



Cat. No.	Color	Description	Sterile	Packaging
PK00909	Blue	Disposable Anti-Static Polypropylene Micro-Sample Transfer Scoop.	YES	100pieces/pack. 60packs/carton

01 SAMPLE HANDLING



MICRO PLASTIC TUBE

PEKYBIO Micro Plastic Tubes, manufactured from high-purity polypropylene, feature a secure screw-cap with integrated sealing ring. Engineered for leak-free performance during centrifugation, storage, and transport, these tubes prevent sample contamination while maintaining chemical resistance. Compatible with ceramic beads for protein and nucleic acid extraction.

- ◆ Conical bottom design for efficient sample recovery.
- ◆ Color-coded caps for sample identification.
- ◆ The attached cap mechanism prevents cap loss.
- ◆ Multiple design options are available.
- ◆ Certified RNase/ DNase free and Non-pyrogen.



PK01215

PK01190

PK00716

PK01236

PK01187

PK00867

ceramic beads



Micro Plastic Tube



Cat. No.	Volume	Color of Tube	Color of Cap	Type of Cap	O-ring	Type of Bottom	Sterile	Packaging
Cap and Tube Pre-assembly								
PK00716	2 mL	Natural	Natural	Screw Cap	Black Silicone	Self Standing	NO	500pieces/bag. 10bags/carton
PK00883	2 mL	Natural	Natural	Screw Cap	Black Silicone	Self Standing	YES	500pieces/bag. 10bags/carton
PK01236	2 mL	Amber	Amber	Screw Cap	Black Silicone	Self Standing	NO	500pieces/bag. 10bags/carton
PK01237	2 mL	Amber	Amber	Screw Cap	Black Silicone	Self Standing	YES	500pieces/bag. 10bags/carton
PK01187	2 mL	Natural	Natural	Screw Cap	Black Silicone	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK01189	2 mL	Natural	Natural	Screw Cap	Black Silicone	Conical Bottom	YES	500pieces/bag. 10bags/carton
Cap and Tube Separation								
PK00715	2 mL	Natural	Natural	Screw Cap	Black Silicone	Self Standing	NO	500pieces/bag. 10bags/carton
PK00925	2 mL	Natural	Natural	Screw Cap	Black Silicone	Self Standing	YES	500pieces/bag. 10bags/carton
PK01085	2 mL	Amber	Amber	Screw Cap	Black Silicone	Self Standing	NO	500pieces/bag. 10bags/carton
PK01235	2 mL	Amber	Amber	Screw Cap	Black Silicone	Self Standing	YES	500pieces/bag. 10bags/carton
PK01186	2 mL	Natural	Natural	Screw Cap	Black Silicone	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK01188	2 mL	Natural	Natural	Screw Cap	Black Silicone	Conical Bottom	YES	500pieces/bag. 10bags/carton
Cap and Tube Attachment								
PK01215	2 mL	Natural	Natural	Screw Loop Cap	Black Silicone	Self Standing	NO	500pieces/bag. 10bags/carton
PK01234	2 mL	Natural	Natural	Screw Loop Cap	Black Silicone	Self Standing	YES	500pieces/bag. 10bags/carton
PK01190	2 mL	Natural	Natural	Screw Loop Cap	Black Silicone	Conical Bottom	NO	500pieces/bag. 10bags/carton
PK01191	2 mL	Natural	Natural	Screw Loop Cap	Black Silicone	Conical Bottom	YES	500pieces/bag. 10bags/carton

Micro Plastic Tube with 1.4mm Ceramic Beads



Cat. No.	Volume	Color of Tube	Color of Cap	Type of Cap	O-ring	Type of Bottom	Sterile	Packaging
PK00867	2 mL	Natural	Natural	Screw Cap	Black Silicone	Self Standing	NO	500pieces/bag. 10bags/carton

◆ All listed tubes are compatible with ceramic beads.

01 SAMPLE HANDLING



CRYOGENIC STORAGE TUBE

PEKYBIO Cryogenic Storage Tubes are manufactured from polypropylene to withstand temperatures to -196°C , range used for sample processing application for small volume sample, excellent for cryo storage or secure sample handling.

- ◆ Multiple working volumes for choice.
- ◆ Free of detectable human DNA, endotoxins, DNase and RNase, Non-cytotoxic.
- ◆ Stored to be used at a temperature between -196°C to 121°C .
- ◆ Easy-to-read graduations.
- ◆ Large white portion for easy writing.
- ◆ Round bottom with self-standing, can single-handed manipulation.



1.8 mL Cardboard Box



4.8 mL Cardboard Box



1.8 mL

4.8 mL

Sterile Cryogenic Storage Tube- 1.8, 4.8 mL



Cat. No.	Volume	Color of Cap	Type of Cap	Color of Tube	Threaded Type	Sterile	Packaging
PK01112	1.8 mL	Blue	Screw Cap	Natural	Externally Threaded	YES	50pieces/bag. 100bags/carton
PK00778	4.8 mL	Blue	Screw Cap	Natural	Externally Threaded	YES	20pieces/bag. 100bags/carton
PK01114	1.8 mL	Blue	Screw Cap	Natural	Externally Threaded	YES	81pieces/box. 10boxes/carton
PK01116	4.8 mL	Blue	Screw Cap	Natural	Externally Threaded	YES	81pieces/box. 10boxes/carton

Non-Sterile Cryogenic Storage Tube- 1.8, 4.8 mL



Cat. No.	Volume	Color of Cap	Type of Cap	Color of Tube	Threaded Type	Sterile	Packaging
PK00540	1.8 mL	Blue	Screw Cap	Natural	Externally Threaded	NO	50pieces/bag. 100bags/carton
PK00541	4.8 mL	Blue	Screw Cap	Natural	Externally Threaded	NO	20pieces/bag. 100bags/carton
PK01113	1.8 mL	Blue	Screw Cap	Natural	Externally Threaded	NO	81pieces/box. 10boxes/carton
PK01115	4.8 mL	Blue	Screw Cap	Natural	Externally Threaded	NO	81pieces/box. 10boxes/carton

Empty Cryogenic Cardboard Box- 81 Hole



Cat. No.	Description	Sterile	Packaging
PK01117	81Hole Cryogenic Cardboard Box for 1.8mL Tube.	NO	5pieces/pack. 10packs/carton
PK01118	81Hole Cryogenic Cardboard Box for 4.8mL Tube.	NO	3pieces/pack. 10packs/carton

01 SAMPLE HANDLING



PCR TUBE & PCR PLATE

PEKYBIO PCR Tubes and Plates are manufactured from premium raw materials and are suitable for a wide range of PCR applications. The product line includes single tubes, tube strips, cap strips, and 96-well PCR plates. All manufactured without chemical additives that could contaminate reagents or samples. PCR Plates are available in multiple formats, including no-skirt and semi-skirt designs. These products are widely used in applications such as genotyping, cloning, mutation detection, sequencing, microarrays, forensics and paternity testing.



HIGH SUITABILITY

**UNIFORM
WALL THICKNESS**



**HIGH-CONTRAST BLACK
PRINTED MARKINGS**

HIGH TRANSPARENCY



8-Strip PCR Tube

- ◆ Compatible with a wide variety of thermal cyclers for parallel PCR workflows.
- ◆ Ultra-thin and consistent wall thickness enable uniform thermal transfer across all tubes.
- ◆ Domed or flat caps ensure a tight seal and effectively prevent sample evaporation.
- ◆ Ideal for multi-sample processing and high-throughput PCR applications.
- ◆ Certified RNase / DNase-free and pyrogen-free.



8-Strip PCR Tube- 0.1, 0.2 mL



Cat. No.	Volume	Description	Type of Cap	Color of Tube	Color of Cap	Sterile	Packaging
PK00187	0.1 mL	8-Strip PCR Tube	Separate Flat Cap	Clear	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton
PK00188	0.1 mL	8-Strip PCR Tube	Separate Flat Cap	White	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton
PK00189	0.1 mL	8-Strip PCR Tube	Separate Domed Cap	Clear	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton
PK00190	0.1 mL	8-Strip PCR Tube	Separate Domed Cap	White	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton
PK00191	0.2 mL	8-Strip PCR Tube	Separate Flat Cap	Clear	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton
PK00192	0.2 mL	8-Strip PCR Tube	Separate Flat Cap	White	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton
PK00193	0.2 mL	8-Strip PCR Tube	Separate Domed Cap	Clear	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton
PK00194	0.2 mL	8-Strip PCR Tube	Separate Domed Cap	White	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton
PK00195	0.1 mL	8-Strip PCR Tube	Attached Flat Cap	Clear	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton
PK00196	0.2 mL	8-Strip PCR Tube	Attached Flat Cap	Clear	Clear	NO	125strips/bag. 1bag/pack. 10packs/carton

Single PCR Tube

- Individual PCR tubes feature evaporation resistance.
- Easy-to-open and close caps prevent sample loss.
- Ultra-thin, uniform tube walls ensure optimal heat transfer and reaction efficiency.
- Sealed in a user-friendly, resealable bag to maintain freshness and minimize contamination.
- Certified free from detectable RNase, DNase, PCR inhibitors and pyrogen.



Single PCR Tube - 0.2, 0.5 mL



Cat. No.	Volume	Description	Type of Cap	Color	Sterile	Packaging
PK00178	0.2 mL	Single PCR Tube	Attached Flat Cap	Clear	NO	1000pieces/bag.10bags/carton
PK00179	0.2 mL	Single PCR Tube	Attached Domed Cap	Clear	NO	1000pieces/bag.10bags/carton
PK00180	0.5 mL	Single PCR Tube	Attached Flat Cap	Clear	NO	1000pieces/bag.10bags/carton

01 SAMPLE HANDLING

PCR Plate



UNIFORM WALL THICKNESS



MULTIPLE SPECS



- ◆ Made of premium polypropylene raw material.
- ◆ Certified free from DNase, RNase and pyrogen.
- ◆ Ultra-thin wall design ensures rapid and uniform heat transfer.
- ◆ Compatible with most thermal cyclers.

PCR Plate- 0.1, 0.2 mL



Cat. No.	Volume	Description	Skirt Type	Color	Deck Type	Sterile	Packaging
PK00197	0.1 mL	96-Well PCR Plate	No Skirt	Frost	Flat	NO	5pieces/bag. 5bags/pack. 6packs/carton
PK00198	0.1 mL	96-Well PCR Plate	No Skirt	White	Flat	NO	5pieces/bag. 5bags/pack. 6packs/carton
PK00199	0.1 mL	96-Well PCR Plate	Semi Skirt	Frost	Non-Flat	NO	5pieces/bag. 5bags/pack. 6packs/carton
PK00200	0.1 mL	96-Well PCR Plate	Semi Skirt	White	Non-Flat	NO	5pieces/bag. 5bags/pack. 6packs/carton
PK00201	0.2 mL	96-Well PCR Plate	No Skirt	Clear	Flat	NO	10pieces/bag. 1bag/pack. 5packs/carton
PK00202	0.2 mL	96-Well PCR Plate	No Skirt	White	Flat	NO	10pieces/bag. 1bag/pack. 5packs/carton
PK00203	0.2 mL	96-Well PCR Plate	Semi Skirt	Clear	Flat	NO	10pieces/bag. 1bag/pack. 5packs/carton
PK00204	0.2 mL	96-Well PCR Plate	Semi Skirt	White	Flat	NO	10pieces/bag. 1bag/pack. 5packs/carton

Sealing Film



EASY TO TEAR



HIGH TRANSPARENCY



- ◆ Compatible with most PCR plates.
- ◆ Non-permeable film, highly transparent.
- ◆ Easy-to-tear design for quick sealing and clean removal.
- ◆ Reliable sealing performance for thermal cycling, excellent optical clarity.

Sealing Film for PCR Plate



Cat. No.	Description	Dimension	Color	Sterile	Packaging
PK00205	PP Sealing Film	141.5 mm x 80 mm	Clear	NO	100pieces/bag. 1bag/pack. 100packs/carton
PK00537	PE Sealing Film	141.5 mm x 80 mm	Clear	NO	100pieces/bag. 1bag/pack. 100packs/carton
PK00539	PET Sealing Film	125mm x 80mm	Clear	NO	100pieces/bag. 1bag/pack. 100packs/carton
PK00538	PP Sealing Paddle for Sealing Film	/	Blue	NO	10pieces/bag

- ◆ PP (PK00205) and PE (PK00537) sealing films are used with sealing paddle; PET sealing film (PK00539) is used in heat sealer.



DEEP WELL PLATE

PEKYBIO Deep Well Plates are designed for compatibility with common laboratory instruments and can be used for applications such as sample collection, compound preparation, combinatorial chemistry, high-throughput screening, nucleic acid purification and bacterial culture.

- ◆ "V"-bottom design permits complete liquid removal.
- ◆ Suitable for long-term sample storage.
- ◆ Meets international SBS standards, compatible with most instruments.
- ◆ Resistant to most chemicals and solvents, including acidic and alkaline solutions.



COMB SET

SBS STANDARD



Deep Well Plate



Cat. No.	Volume	Well	Well Type	Bottom Type	Sterile	Packaging
PK00158	0.36 mL	96-Well	Round	V Bottom	NO	10pieces/bag. 1bag/pack. 10packs/carton
PK00159	0.40 mL	96-Well	Round	U Bottom	NO	10pieces/bag. 1bag/pack. 10packs/carton
PK00160	0.50 mL	96-Well	Square	V Bottom	NO	10pieces/bag. 1bag/pack. 10packs/carton
PK00161	1.00 mL	96-Well	Round	U Bottom	NO	5pieces/bag. 2bags/pack. 5packs/carton
PK00162	1.00 mL	96-Well	Square	U Bottom	NO	10pieces/bag. 2bags/pack. 5packs/carton
PK00163	1.60 mL	96-Well	Square	U Bottom	NO	5pieces/bag. 2bags/pack. 5packs/carton
PK00164	2.00 mL	96-Well	Round	U Bottom	NO	5pieces/bag. 2bags/pack. 5packs/carton

01 SAMPLE HANDLING

Deep Well Plate



Cat. No.	Volume	Well	Well Type	Bottom Type	Sterile	Packaging
PK00165	2.20 mL	96-Well	Square	V Bottom	NO	5pieces/bag. 2bags/pack. 5packs/carton
PK00166	96-Tip Comb	96-Well	Compatible with PK00165		NO	2pieces/bag. 25bags/carton
PK00167	2.20 mL	96-Well	Square	U Bottom	NO	5pieces/bag. 2bags/pack. 5packs/carton
PK00536	8-Strip Tip Comb	8-Well	Compatible with PK00167		NO	2pieces/bag. 100bags/carton
PK00168	2.30 mL	96-Well	Square	V Bottom	NO	5pieces/bag. 2bags/pack. 5packs/carton
PK00169	3.50 mL	48-Well	Round	U Bottom	NO	5pieces/bag. 2bags/pack. 5packs/carton
PK00170	4.60 mL	48-Well	Square	U Bottom	NO	5pieces/bag. 2bags/pack. 5packs/carton
PK00171	10.00 mL	24-Well	Square	V Bottom	NO	5pieces/bag. 2bags/pack. 5packs/carton

SILICONE SEALING MAT

Silicone Sealing Mat



Cat. No.	Well	Pierceable	Thinned-Web	Compatible	Sterile	Packaging
PK00172	96-Well	NO	/	Fit 96-Square-Well	NO	10pieces/bag. 5bags/carton
PK00173	96-Well	YES	/	Fit 96-Square-Well	NO	10pieces/bag. 5bags/carton
PK00174	96-Well	YES	Straight	Fit 2mL 96-Round-Well	NO	10pieces/bag. 5bags/carton
PK00394	96-Well	YES	Cross	Fit 2mL 96-Round-Well	NO	10pieces/bag. 5bags/carton
PK00175	96-Well	YES	Cross	Fit 1mL 96-Round-Well	NO	10pieces/bag. 5bags/carton
PK00176	48-Well	YES	/	Fit 48-Square-Well	NO	10pieces/bag. 5bags/carton



ALUMINUM FOIL SEALING FILM

Aluminum Foil Sealing Film



Cat. No.	Type of Surface	Compatible	Sterile	Packaging
PK00177	Smooth Surface	Self Adhesive	NO	100pieces/bag. 1bag/pack. 100packs/carton
PK00562	Smooth Surface	Used In Heat Sealer	NO	100pieces/bag. 1bag/pack. 100packs/carton
PK00563	Embossed Surface	Used In Heat Sealer	NO	100pieces/bag. 1bag/pack. 100packs/carton



PK00562



PK00563





02 LIQUID HANDLING

Serological Pipette/ [26-27](#)

Aspirating Pipette/ [28](#)

Conductive Automation Tip/ [29](#)

Universal Pipette Tip/ [30-36](#)

LCTS Pipette Tip/ [37-40](#)

Pasteur Pipette/ [41-42](#)

Exact Volume Transfer Pipette/ [42](#)

02 LIQUID HANDLING



SEROLOGICAL PIPETTE

PEKYBIO Serological Pipettes are frequently used in laboratories to transfer liquids. These pipettes have graduations on the side to help measure liquids to be dispensed or aspirated. They are manufactured from optically clear polystyrene in a class 100,000 cleanroom. All pipettes and individual packaging are color-coded for easy volume identification. Available in 1mL, 2mL, 5mL, 10mL, 25mL and 50mL. Peel-to-open plastic/paper wrappers are the ideal option for eliminating particulate.



INDIVIDUAL PACKAGING



ACCURATE DROP




CLEAR GRADUATIONS



OBVIOUS COLOR CODES

- ◆ Clear and precise graduations for easy volume determination.
- ◆ Compatible with most pipettes.
- ◆ Certified free from DNase, RNase and Non-pyrogen.
- ◆ Sterilized by irradiation, SAL=10⁻⁶.
- ◆ Individually wrapped pipettes for reliable sterility.
- ◆ Both the pipettes and individual packaging are color-coded for easy volume identification.



1 mL

2 mL

5 mL

10 mL

25 mL

50 mL

Plastic/Paper Individual Packaging - 1, 2, 5, 10, 25, 50 mL



Cat. No.	Volume	Color Code	Sterile	Packaging
PK00086	1 mL	Yellow	YES	100pieces/bag. 3bags/pack. 6packs/carton
PK00087	2 mL	Green	YES	100pieces/bag. 3bags/pack. 6packs/carton
PK00088	5 mL	Blue	YES	50pieces/bag. 4bags/pack. 6packs/carton
PK00089	10 mL	Orange	YES	50pieces/bag. 4bags/pack. 6packs/carton
PK00090	25 mL	Red	YES	25pieces/bag. 4bags/pack. 6packs/carton
PK00091	50 mL	Purple	YES	25pieces/bag. 3bags/pack. 6packs/carton

Plastic Individual Packaging - 1, 2, 5, 10, 25, 50 mL



Cat. No.	Volume	Color Code	Sterile	Packaging
PK01130	1 mL	Yellow	YES	100pieces/bag. 3bags/pack. 6packs/carton
PK01131	2 mL	Green	YES	100pieces/bag. 3bags/pack. 6packs/carton
PK01132	5 mL	Blue	YES	50pieces/bag. 4bags/pack. 6packs/carton
PK01133	10 mL	Orange	YES	50pieces/bag. 4bags/pack. 6packs/carton
PK01134	25 mL	Red	YES	25pieces/bag. 4bags/pack. 6packs/carton
PK01135	50 mL	Purple	YES	25pieces/bag. 3bags/pack. 6packs/carton



Plastic Bag of 25pcs Packaging - 1, 2, 5, 10, 25, 50 mL



Cat. No.	Volume	Color Code	Sterile	Packaging
PK00659	1 mL	Yellow	YES	25pieces/bag. 20bags/pack. 3packs/carton
PK00660	2 mL	Green	YES	25pieces/bag. 16bags/pack. 3packs/carton
PK00661	5 mL	Blue	YES	25pieces/bag. 10bags/pack. 6packs/carton
PK00662	10 mL	Orange	YES	25pieces/bag. 8bags/pack. 6packs/carton
PK00663	25 mL	Red	YES	25pieces/bag. 5bags/pack. 6packs/carton
PK00664	50 mL	Purple	YES	25pieces/bag. 4bags/pack. 6packs/carton

ASPIRATING PIPETTE

PEKYBIO Aspirating Pipettes are completely transparent and designed for the transfer and mixing of liquids. It is a safe alternative to glass pipettes. It is also used in a wide variety of applications such as vacuum aspiration or pipettor-based aspiration procedures.

- ◆ Ideal for tissue culture applications and liquid transfer.
- ◆ Sterile, DNase / RNase-free, and Non-pyrogenic.
- ◆ Individually wrapped in paper-plastic packaging.



Aspirating Pipette - 1, 2, 5, 10, 25, 50 mL



Cat. No.	Volume	Packing Status	Sterile	Packaging	
PK00707	1 mL	Plastic/Paper Individual Packaging.	YES	1800 pcs	100pieces/bag. 3bags/pack. 6packs/carton
PK00708	2 mL	Plastic/Paper Individual Packaging.	YES	1800 pcs	100pieces/bag. 3bags/pack. 6packs/carton
PK00709	5 mL	Plastic/Paper Individual Packaging.	YES	1200 pcs	50pieces/bag. 4bags/pack. 6packs/carton
PK00710	10 mL	Plastic/Paper Individual Packaging.	YES	1200 pcs	50pieces/bag. 4bags/pack. 6packs/carton
PK00711	25 mL	Plastic/Paper Individual Packaging.	YES	600 pcs	25pieces/bag. 4bags/pack. 6packs/carton
PK00712	50 mL	Plastic/Paper Individual Packaging.	YES	450 pcs	25pieces/bag. 3bags/pack. 6packs/carton

CONDUCTIVE AUTOMATION TIP

Automated pipetting processes place high demands on the design and raw materials of pipette tips. Filter pipette tip helps to reduce the risk of aerosols forming in the experiments. It also protects the pipette shafts from contamination and reduces the risk of cross contamination. The pipette tip is often used in contamination sensitive applications such as forensics and clinical diagnostics.



INDIVIDUAL RACK



BLACK AND CONDUCTIVE



HIGHLY COMPATIBLE



MULTIPLE VOLUMES

- ◆ Conductive automation tips are subject to tighter tolerances than tips designed for handheld pipetting.
- ◆ 100% automated assembly to ensure consistent performance and quality.
- ◆ Compatible with Hamilton and Tecan liquid handling platforms.

Rack Packaging ◆ 96 tips per rack. ◆ 8 x 12 format for use with multi-channel pipettors.

Conductive Automation Tip- 50, 200, 300, 1000 uL



Cat. No.	Volume	Filter	Color	Sterile	Packaging
Compatible with Hamilton					
PKH0015	50 uL	NO	Black	NO	96pieces/rack. 24racks/carton
PKH0017	50 uL	NO	Black	YES	96pieces/rack. 24racks/carton
PKH0016	50 uL	YES	Black	NO	96pieces/rack. 24racks/carton
PKH0018	50 uL	YES	Black	YES	96pieces/rack. 24racks/carton
PKH0019	300 uL	NO	Black	NO	96pieces/rack. 24racks/carton
PKH0021	300 uL	NO	Black	YES	96pieces/rack. 24racks/carton
PKH0020	300 uL	YES	Black	NO	96pieces/rack. 24racks/carton
PKH0022	300 uL	YES	Black	YES	96pieces/rack. 24racks/carton
PKH0023	1000 uL	NO	Black	NO	96pieces/rack. 24racks/carton
PKH0025	1000 uL	NO	Black	YES	96pieces/rack. 24racks/carton
PKH0024	1000 uL	YES	Black	NO	96pieces/rack. 24racks/carton
PKH0026	1000 uL	YES	Black	YES	96pieces/rack. 24racks/carton
Compatible with Tecan					
PKT0027	50 uL	NO	Black	NO	96pieces/rack. 24racks/carton
PKT0029	50 uL	NO	Black	YES	96pieces/rack. 24racks/carton
PKT0028	50 uL	YES	Black	NO	96pieces/rack. 24racks/carton
PKT0030	50 uL	YES	Black	YES	96pieces/rack. 24racks/carton
PKT0031	200 uL	NO	Black	NO	96pieces/rack. 24racks/carton
PKT0033	200 uL	NO	Black	YES	96pieces/rack. 24racks/carton
PKT0032	200 uL	YES	Black	NO	96pieces/rack. 24racks/carton
PKT0034	200 uL	YES	Black	YES	96pieces/rack. 24racks/carton
PKT0035	1000 uL	NO	Black	NO	96pieces/rack. 24racks/carton
PKT0037	1000 uL	NO	Black	YES	96pieces/rack. 24racks/carton
PKT0036	1000 uL	YES	Black	NO	96pieces/rack. 24racks/carton
PKT0038	1000 uL	YES	Black	YES	96pieces/rack. 24racks/carton

02 LIQUID HANDLING



UNIVERSAL PIPETTE TIP

PEKYBIO Pipette Tips are the most widely used consumables in laboratories, including chemical and biological labs. PEKYBIO high precision pipette tip is guaranteed by precise molding and top-ranking injection machines. Filter pipette tip helps to reduce the risk of aerosols forming in the experiments. They also protect the pipette shafts from contamination and reduce the risk of cross contamination. They are often used in contamination sensitive applications such as forensics and clinical diagnostics.

- ◆ Made from medical-grade virgin polypropylene (PP).
- ◆ Certified RNase free, DNase free, human DNA-free and Non-Pyrogen.
- ◆ Reliable protection from liquids.
- ◆ Compatible with most brands of pipettors.

Rack Packaging

- ◆ 96 tips per rack.
- ◆ 8x12 format for use with single and multi-channel pipettors.

Bag Packaging

- ◆ Economical choice.
- ◆ Each bag features a tamper-evident seal and is resealable to keep out dust and particles.

Refill Packaging

- ◆ Ideal for space-saving and convenient storage.
- ◆ Each package is designed for easy access and efficient handling, ensuring a neat and organized lab environment.

Rack Packaging

Universal Pipette Tip- Rack Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK00001	10 uL	NO	Clear	31 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00599	10 uL	NO	Clear	45 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00004	200 uL	NO	Clear	50.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton

Universal Pipette Tip- Rack Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK00344	200 uL	NO	Yellow	50.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00005	300 uL	NO	Clear	59 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00605	300 uL	NO	Clear	89 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00006	1000 uL	NO	Clear	87 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00006-78	1000 uL	NO	Clear	78 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00395	1000 uL	NO	Clear	78 mm	YES	100pieces/rack. 10racks/pack. 5packs/carton
	PK00345	1000 uL	NO	Blue	87 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00345-78	1000 uL	NO	Blue	78 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00007	1250 uL	NO	Clear	100 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00417	1250 uL	NO	Blue	100 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00759	5 mL	NO	Clear	130 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK00719	5 mL	NO	Clear	150 mm	YES	40pieces/rack. 8racks/pack. 2packs/carton
	PK00723	10 mL	NO	Clear	161 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK00763	10 mL	NO	Clear	165 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK00008	10 uL	YES	Clear	31 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00009	10 uL	YES	Clear	45 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00010	20 uL	YES	Clear	50.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00404	50 uL	YES	Clear	50.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00011	100 uL	YES	Clear	50.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00012	200 uL	YES	Clear	59 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00390	200 uL	YES	Yellow	59 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00571	300 uL	YES	Clear	59 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00549	300 uL	YES	Clear	89 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00013	1000 uL	YES	Clear	87 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00013-78	1000 uL	YES	Clear	78 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00396	1000 uL	YES	Clear	78 mm	YES	100pieces/rack. 10racks/pack. 5packs/carton
	PK00393	1000 uL	YES	Blue	87 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00393-78	1000 uL	YES	Blue	78 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00014	1250 uL	YES	Clear	100 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00418	1250 uL	YES	Blue	100 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00760	5 mL	YES	Clear	130 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK00720	5 mL	YES	Clear	150 mm	YES	40pieces/rack. 8racks/pack. 2packs/carton
	PK00724	10 mL	YES	Clear	161 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK00764	10 mL	YES	Clear	165 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton


Bag Packaging

Universal Pipette Tip- Bag Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK00039	10 uL	NO	Clear	31 mm	NO	1000pieces/bag. 20bags/carton
	PK00557	10 uL	NO	Clear	45 mm	NO	1000pieces/bag. 20bags/carton
	PK00060	200 uL	NO	Clear	50.5 mm	NO	1000pieces/bag. 20bags/carton
	PK00388	200 uL	NO	Yellow	50.5 mm	NO	1000pieces/bag. 20bags/carton
	PK00066	300 uL	NO	Clear	59 mm	NO	1000pieces/bag. 20bags/carton
	PK00443	300 uL	NO	Clear	89 mm	NO	1000pieces/bag. 10bags/carton
	PK00069	1000 uL	NO	Clear	87 mm	NO	1000pieces/bag. 5bags/carton
	PK00069-78	1000 uL	NO	Clear	78 mm	NO	1000pieces/bag. 5bags/carton
	PK00389	1000 uL	NO	Blue	87 mm	NO	1000pieces/bag. 5bags/carton
	PK00389-78	1000 uL	NO	Blue	78 mm	NO	1000pieces/bag. 5bags/carton
	PK00075	1250 uL	NO	Clear	100 mm	NO	1000pieces/bag. 5bags/carton
	PK00419	1250 uL	NO	Blue	100 mm	NO	1000pieces/bag. 5bags/carton
	PK00078	5 mL	NO	Clear	130 mm	NO	100pieces/bag. 10bags/carton
	PK00375	5 mL	NO	Clear	150 mm	NO	100pieces/bag. 10bags/carton
	PK00082	10 mL	NO	Clear	161 mm	NO	100pieces/bag. 10bags/carton
	PK00379	10 mL	NO	Clear	165 mm	NO	100pieces/bag. 10bags/carton
	PK00040	10 uL	YES	Clear	31 mm	NO	1000pieces/bag. 20bags/carton
	PK00041	10 uL	YES	Clear	45 mm	NO	1000pieces/bag. 20bags/carton
	PK00049	20 uL	YES	Clear	50.5 mm	NO	1000pieces/bag. 20bags/carton
	PK00405	50 uL	YES	Clear	50.5 mm	NO	1000pieces/bag. 20bags/carton
	PK00055	100 uL	YES	Clear	50.5 mm	NO	1000pieces/bag. 20bags/carton
	PK00061	200 uL	YES	Clear	59 mm	NO	1000pieces/bag. 20bags/carton
	PK00615	200 uL	YES	Yellow	59 mm	NO	1000pieces/bag. 20bags/carton
	PK00602	300 uL	YES	Clear	59 mm	NO	1000pieces/bag. 20bags/carton
	PK00601	300 uL	YES	Clear	89 mm	NO	1000pieces/bag. 10bags/carton
	PK00070	1000 uL	YES	Clear	87 mm	NO	1000pieces/bag. 5bags/carton
	PK00070-78	1000 uL	YES	Clear	78 mm	NO	1000pieces/bag. 5bags/carton
	PK00416	1000 uL	YES	Blue	87 mm	NO	1000pieces/bag. 5bags/carton
	PK00416-78	1000 uL	YES	Blue	78 mm	NO	1000pieces/bag. 5bags/carton
	PK00398	1250 uL	YES	Clear	100 mm	NO	1000pieces/bag. 5bags/carton
	PK00420	1250 uL	YES	Blue	100 mm	NO	1000pieces/bag. 5bags/carton
	PK00080	5 mL	YES	Clear	130 mm	NO	100pieces/bag. 10bags/carton
	PK00377	5 mL	YES	Clear	150 mm	NO	100pieces/bag. 10bags/carton
	PK00381	10 mL	YES	Clear	161 mm	NO	100pieces/bag. 10bags/carton
	PK00084	10 mL	YES	Clear	165 mm	NO	100pieces/bag. 10bags/carton

Refill Packaging



Universal Pipette Tip- Refill Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK01196	10 uL	NO	Clear	45 mm	YES	960pieces/pack. 10packs/carton
	PK01197	200 uL	NO	Clear	50.5 mm	YES	960pieces/pack. 10packs/carton
	PK01198	200 uL	NO	Yellow	50.5 mm	YES	960pieces/pack. 10packs/carton
	PK01199	300 uL	NO	Clear	59 mm	YES	960pieces/pack. 10packs/carton
	PK01200	1000 uL	NO	Clear	87 mm	YES	960pieces/pack. 5packs/carton
	PK01201	1000 uL	NO	Blue	87 mm	YES	960pieces/pack. 5packs/carton
	PK01202	1250 uL	NO	Clear	100 mm	YES	960pieces/pack. 5packs/carton
	PK01203	1250 uL	NO	Blue	100 mm	YES	960pieces/pack. 5packs/carton

Low Retention - Rack Packaging

Universal Pipette Tip- Low Retention- Rack Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK00045	10 uL	NO	Clear	31 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00600	10 uL	NO	Clear	45 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00064	200 uL	NO	Clear	50.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00068	300 uL	NO	Clear	59 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00606	300 uL	NO	Clear	89 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00073	1000 uL	NO	Clear	87 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00073-78	1000 uL	NO	Clear	78 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01091	1000 uL	NO	Clear	78 mm	YES	100pieces/rack. 10racks/pack. 5packs/carton
	PK00077	1250 uL	NO	Clear	100 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00761	5 mL	NO	Clear	130 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK00721	5 mL	NO	Clear	150 mm	YES	40pieces/rack. 8racks/pack. 2packs/carton
	PK00725	10 mL	NO	Clear	161 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK00765	10 mL	NO	Clear	165 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton

02 LIQUID HANDLING

Universal Pipette Tip- Low Retention- Rack Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK00046	10uL	YES	Clear	31 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00047	10 uL	YES	Clear	45 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00053	20 uL	YES	Clear	50.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00407	50 uL	YES	Clear	50.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00059	100 uL	YES	Clear	50.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00065	200 uL	YES	Clear	59 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00607	300 uL	YES	Clear	59 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00617	300 uL	YES	Clear	89 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00074	1000 uL	YES	Clear	87 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00074-78	1000 uL	YES	Clear	78 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01092	1000 uL	YES	Clear	78 mm	YES	100pieces/rack. 10racks/pack. 5packs/carton
	PK00397	1250 uL	YES	Clear	100 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK00762	5 mL	YES	Clear	130 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK00722	5 mL	YES	Clear	150 mm	YES	40pieces/rack. 8racks/pack. 2packs/carton
	PK00726	10 mL	YES	Clear	161 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK00766	10 mL	YES	Clear	165 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton

Low Retention - Bag Packaging

Universal Pipette Tip- Low Retention- Bag Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK00042	10 uL	NO	Clear	31 mm	NO	1000pieces/bag. 20bags/carton
	PK00598	10 uL	NO	Clear	45 mm	NO	1000pieces/bag. 20bags/carton
	PK00062	200 uL	NO	Clear	50.5 mm	NO	1000pieces/bag. 20bags/carton
	PK00067	300 uL	NO	Clear	59 mm	NO	1000pieces/bag. 20bags/carton
	PK00603	300 uL	NO	Clear	89 mm	NO	1000pieces/bag. 10bags/carton
	PK00071	1000 uL	NO	Clear	87 mm	NO	1000pieces/bag. 5bags/carton
	PK00071-78	1000 uL	NO	Clear	78 mm	NO	1000pieces/bag. 5bags/carton
	PK00076	1250 uL	NO	Clear	100 mm	NO	1000pieces/bag. 5bags/carton
	PK00079	5 mL	NO	Clear	130 mm	NO	100pieces/bag. 10bags/carton
	PK00376	5 mL	NO	Clear	150 mm	NO	100pieces/bag. 10bags/carton
	PK00083	10 mL	NO	Clear	161 mm	NO	100pieces/bag. 10bags/carton
	PK00380	10 mL	NO	Clear	165 mm	NO	100pieces/bag. 10bags/carton
	PK00043	10 uL	YES	Clear	31 mm	NO	1000pieces/bag. 20bags/carton
	PK00044	10 uL	YES	Clear	45 mm	NO	1000pieces/bag. 20bags/carton
	PK00051	20 uL	YES	Clear	50.5 mm	NO	1000pieces/bag. 20bags/carton
	PK00406	50 uL	YES	Clear	50.5 mm	NO	1000pieces/bag. 20bags/carton

Low Retention - Bag Packaging

Universal Pipette Tip- Low Retention- Bag Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK00057	100 uL	YES	Clear	50.5 mm	NO	1000pieces/bag. 20bags/carton
	PK00063	200 uL	YES	Clear	59 mm	NO	1000pieces/bag. 20bags/carton
	PK00604	300 uL	YES	Clear	59 mm	NO	1000pieces/bag. 20bags/carton
	PK00616	300 uL	YES	Clear	89 mm	NO	1000pieces/bag. 10bags/carton
	PK00072	1000 uL	YES	Clear	87 mm	NO	1000pieces/bag. 5bags/carton
	PK00072-78	1000 uL	YES	Clear	78 mm	NO	1000pieces/bag. 5bags/carton
	PK00399	1250 uL	YES	Clear	100 mm	NO	1000pieces/bag. 5bags/carton
	PK00081	5 mL	YES	Clear	130 mm	NO	100pieces/bag. 10bags/carton
	PK00378	5 mL	YES	Clear	150 mm	NO	100pieces/bag. 10bags/carton
	PK00085	10 mL	YES	Clear	161 mm	NO	100pieces/bag. 10bags/carton
	PK00382	10 mL	YES	Clear	165 mm	NO	100pieces/bag. 10bags/carton


Low Retention - Refill Packaging

Universal Pipette Tip- Refill Packaging




Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK01204	10 uL	NO	Clear	45 mm	YES	960pieces/pack. 10packs/carton
	PK01205	200 uL	NO	Clear	50.5 mm	YES	960pieces/pack. 10packs/carton
	PK01206	300 uL	NO	Clear	59 mm	YES	960pieces/pack. 10packs/carton
	PK01207	1000 uL	NO	Clear	87 mm	YES	960pieces/pack. 5packs/carton
	PK01208	1250 uL	NO	Clear	100 mm	YES	960pieces/pack. 5packs/carton

Length for Reference Only



02 LIQUID HANDLING





LCTS PIPETTE TIP

PEKYBIO Light Contact Tip Solution (LCTS) Pipette Tips ensure exceptional accuracy, reliability, and effortless tip mounting and ejection, making them ideal for precise laboratory applications. Available in multiple volumes, they are perfect for molecular biology, cell culture and other critical liquid handling tasks.

- ◆ Excellent sealing and easy ejection with minimal force.
- ◆ Minimizes sample loss, maximizing precision and efficiency in experiments.
- ◆ Certified free of RNase, DNase, Human DNA, ATP, and endotoxins for reliable results.
- ◆ Clear graduations for easy volume matching and accurate measurements.
- ◆ Made from virgin polypropylene, providing outstanding clarity and chemical resistance.

Rack Packaging

LCTS Pipette Tip- Rack Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK01008	20 uL	NO	Clear	49.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01009	200 uL	NO	Clear	50 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01010	300 uL	NO	Clear	57.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01011	1000 uL	NO	Clear	86 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01012	1250 uL	NO	Clear	86 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01013	5 mL	NO	Clear	130 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK01014	10 mL	NO	Clear	165 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton

02 LIQUID HANDLING

LCTS Pipette Tip- Rack Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK01015	20 uL	YES	Clear	49.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01016	200 uL	YES	Clear	50 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01017	300 uL	YES	Clear	57.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01018	1000 uL	YES	Clear	86 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01019	1250 uL	YES	Clear	86 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01020	5 mL	YES	Clear	130 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK01021	10 mL	YES	Clear	165 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton

Bag Packaging

LCTS Pipette Tip- Bag Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK01022	20 uL	NO	Clear	49.5 mm	NO	1000pieces/bag. 10bags/carton
	PK01023	200 uL	NO	Clear	50 mm	NO	1000pieces/bag. 10bags/carton
	PK01024	300 uL	NO	Clear	57.5 mm	NO	1000pieces/bag. 10bags/carton
	PK01025	1000 uL	NO	Clear	86 mm	NO	1000pieces/bag. 10bags/carton
	PK01026	1250 uL	NO	Clear	86 mm	NO	1000pieces/bag. 10bags/carton
	PK01027	5 mL	NO	Clear	130 mm	NO	100pieces/bag. 10bags/carton
	PK01028	10 mL	NO	Clear	165 mm	NO	100pieces/bag. 10bags/carton
	PK01029	20 uL	YES	Clear	49.5 mm	NO	1000pieces/bag. 10bags/carton
	PK01030	200 uL	YES	Clear	50 mm	NO	1000pieces/bag. 10bags/carton
	PK01031	300 uL	YES	Clear	57.5 mm	NO	1000pieces/bag. 10bags/carton
	PK01032	1000 uL	YES	Clear	86 mm	NO	1000pieces/bag. 10bags/carton
	PK01033	1250 uL	YES	Clear	86 mm	NO	1000pieces/bag. 10bags/carton
	PK01034	5 mL	YES	Clear	130 mm	NO	100pieces/bag. 10bags/carton
	PK01035	10 mL	YES	Clear	165 mm	NO	100pieces/bag. 10bags/carton

Low Retention - Rack Packaging

LCTS Pipette Tip- Low Retention- Rack Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK01036	20 uL	NO	Clear	49.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01037	200 uL	NO	Clear	50 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01038	300 uL	NO	Clear	57.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01039	1000 uL	NO	Clear	86 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01040	1250 uL	NO	Clear	86 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01041	5 mL	NO	Clear	130 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK01042	10 mL	NO	Clear	165 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK01043	20 uL	YES	Clear	49.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01044	200 uL	YES	Clear	50 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01045	300 uL	YES	Clear	57.5 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01046	1000 uL	YES	Clear	86 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01047	1250 uL	YES	Clear	86 mm	YES	96pieces/rack. 10racks/pack. 5packs/carton
	PK01048	5 mL	YES	Clear	130 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton
	PK01049	10 mL	YES	Clear	165 mm	YES	24pieces/rack. 8racks/pack. 2packs/carton

Low Retention - Bag Packaging

LCTS Pipette Tip- Low Retention- Bag Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK01050	20 uL	NO	Clear	49.5 mm	NO	1000pieces/bag. 10bags/carton
	PK01051	200 uL	NO	Clear	50 mm	NO	1000pieces/bag. 10bags/carton
	PK01052	300 uL	NO	Clear	57.5 mm	NO	1000pieces/bag. 10bags/carton
	PK01053	1000 uL	NO	Clear	86 mm	NO	1000pieces/bag. 10bags/carton
	PK01054	1250 uL	NO	Clear	86 mm	NO	1000pieces/bag. 10bags/carton
	PK01055	5 mL	NO	Clear	130 mm	NO	100pieces/bag. 10bags/carton
	PK01056	10 mL	NO	Clear	165 mm	NO	100pieces/bag. 10bags/carton


02 LIQUID HANDLING

LCTS Pipette Tip- Low Retention- Bag Packaging



Photo	Cat. No.	Volume	Filter	Color	Length	Sterile	Packaging
	PK01057	20 uL	YES	Clear	49.5 mm	NO	1000pieces/bag. 10bags/carton
	PK01058	200 uL	YES	Clear	50 mm	NO	1000pieces/bag. 10bags/carton
	PK01059	300 uL	YES	Clear	57.5 mm	NO	1000pieces/bag. 10bags/carton
	PK01060	1000 uL	YES	Clear	86 mm	NO	1000pieces/bag. 10bags/carton
	PK01061	1250 uL	YES	Clear	86 mm	NO	1000pieces/bag. 10bags/carton
	PK01062	5 mL	YES	Clear	130 mm	NO	100pieces/bag. 10bags/carton
	PK01063	10 mL	YES	Clear	165 mm	NO	100pieces/bag. 10bags/carton

Length for Reference Only





PASTEUR PIPETTE

PEKYBIO Pasteur Pipettes are essential daily tools in the laboratory, designed as a safe and convenient alternative to traditional glass pipettes. Featuring an integrated bulb and single-use design, they help minimize the risk of cross-contamination during routine sample handling. Ideal for general lab applications—including transferring supernatants, adding reagents and collecting samples for analysis.



DISPENSER BOX & STACKABLE DESIGN

INDIVIDUAL PACKAGING




CLEAR GRADUATIONS

SUPERIOR RESILIENCE



- ◆ Non-toxic and chemically inert to biological fluids and most acids.
- ◆ Ideally suited for the rapid transfer of liquids.
- ◆ Non-sterile Pasteur pipettes are packed in a unique dispenser box that is compact and ideal for bench-top use.

Length can be customized.



02 LIQUID HANDLING

Pasteur Pipette- 1, 2, 3, 5, 10 mL



Cat. No.	Volume	Length	Graduation	Color	Sterile	Packaging
Bulk Dispenser Box						
PK00092	1 mL	150 mm	Double Scale	Clear	NO	500pieces/pack. 10packs/carton
PK00093	2 mL	155 mm	Single Scale	Clear	NO	500pieces/pack. 10packs/carton
PK00094	3 mL	155 mm	Single Scale	Clear	NO	500pieces/pack. 10packs/carton
PK00111	3 mL	160 mm	Double Scale	Clear	NO	500pieces/pack. 10packs/carton
Plastic/Paper Individual Packaging						
PK00095	1 mL	150 mm	Double Scale	Clear	YES	500pieces/pack. 4packs/carton
PK00096	2 mL	155 mm	Single Scale	Clear	YES	500pieces/pack. 4packs/carton
PK00097	3 mL	155 mm	Single Scale	Clear	YES	500pieces/pack. 4packs/carton
Plastic Individual Packaging						
PK00696	1 mL	150 mm	Double Scale	Clear	YES	200pieces/bag. 20bags/carton
PK00697	2 mL	155 mm	Single Scale	Clear	YES	200pieces/bag. 20bags/carton
PK00099	3 mL	155 mm	Single Scale	Clear	YES	500pieces/pack. 4packs/carton
PK00698	3 mL	155 mm	Single Scale	Clear	YES	200pieces/bag. 20bags/carton
PK00699	5 mL	210 mm	Single Scale	Clear	YES	100pieces/bag. 30bags/carton
PK01119	10 mL	300 mm	Single Scale	Clear	YES	100pieces/bag. 15bags/carton
Bulk Plastic Bag						
PK00553	1 mL	150 mm	Double Scale	Clear	NO	500pieces/bag. 20bags/carton
PK00554	3 mL	155 mm	Single Scale	Clear	NO	500pieces/bag. 20bags/carton
PK00555	5 mL	210 mm	Single Scale	Clear	NO	100pieces/bag. 50bags/carton
PK00119	10 mL	300 mm	Single Scale	Clear	NO	100pieces/bag. 25bags/carton


EXACT VOLUME TRANSFER PIPETTE

PEKYBIO Dual-Bulb Exact Volume Transfer Pipettes are unbreakable, as they are manufactured from soft plastic which can be recycled when used with non-hazardous material.

Exact volume pipette, also known as Dual-Bulb pasteur pipette is manufactured from low density polyethylene (LDPE) and designed to dispense a known volume of liquid.

This pipette perfect for untrained users as part of a self-test diagnostic kit and it's a quick and easy alternative to pipette and tip.

- Fast, measured transfers.
- Easy to use, economical and accurate.



60 uL PK00120

100 uL PK00121

120 uL PK00122

730 uL PK01109

Exact Volume Transfer Pipette- 60, 100, 120 ,730uL



Cat. No.	Volume	Length	Graduation	Color	Sterile	Packaging From	Packaging
PK00120	60 uL	63 mm	Non-Graduated	Clear	NO	Bulk	500pieces/bag. 40bags/carton
PK00121	100 uL	85 mm	Non-Graduated	Clear	NO	Bulk	500pieces/bag. 40bags/carton
PK00122	120 uL	92 mm	Non-Graduated	Clear	NO	Bulk	500pieces/bag. 40bags/carton
PK01109	730 uL	137 mm	Non-Graduated	Clear	NO	Bulk	100pieces/bag. 100bags/carton



03 CELL CULTURE

Cell Culture Plate- 44-45

Cell Culture Flask- 46-47

Cell Culture Dish- 48-49

Petri Dish- 49-50

Erlenmeyer Flask- 51-52

03 CELL CULTURE



CELL CULTURE PLATE

PEKYBIO Cell Culture Plates are indispensable tools in cellular research, widely used for culturing and observing cells in laboratory and clinical environments. Their flat, smooth surface is designed to optimize cell attachment and growth, making them ideal for cell-based assays and experiments. The uniform hole positions ensure consistent cell distribution, providing reliable and reproducible results. They are made from high-quality polystyrene, ensuring excellent clarity, chemical resistance and reliable performance throughout long-term cell culture applications.



UNIFORM HOLE POSITION

**100%
VIRGIN POLYSTYRENE**



INDIVIDUAL PACKAGING


EASY TO TEAR




- ◆ High clarity, 100 % virgin polystyrene.
- ◆ "TC" refers to vacuum plasma tissue culture surface treatment.
- ◆ Available in 6, 12, 24, 48 and 96 well plates.
- ◆ Sterilized by irradiation, SAL=10⁻⁶.

Cell Culture Plate Sizes


- PEKYBIO plate for cell culture is available in multiple sizes, designed to suit various cell culture




6-Well




12-Well



24-Well




48-Well




96-Well

96 Well Plates are available in flat and U bottoms



Flat-Bottom



U-Bottom

Cell Culture Plate



Cat. No.	Well	Type of Bottom	Color	Surface Treatment	Sterile	Packaging	
PK00254	6-Well	Flat Bottom	Clear	TC Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00255	12-Well	Flat Bottom	Clear	TC Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00256	24-Well	Flat Bottom	Clear	TC Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00257	48-Well	Flat Bottom	Clear	TC Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00258	96-Well	Flat Bottom	Clear	TC Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00259	96-Well	U Bottom	Clear	TC Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00427	6-Well	Flat Bottom	Clear	Non-Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00428	12-Well	Flat Bottom	Clear	Non-Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00429	24-Well	Flat Bottom	Clear	Non-Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00430	48-Well	Flat Bottom	Clear	Non-Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00431	96-Well	Flat Bottom	Clear	Non-Treated	YES	Individual Packaging	1piece/pack. 50packs/carton
PK00432	96-Well	U Bottom	Clear	Non-Treated	YES	Individual Packaging	1piece/pack. 50packs/carton

03 CELL CULTURE



CELL CULTURE FLASK

PEKYBIO Cell Culture Flasks are essential for large-scale cell growth, commonly used in research, biotechnology and pharmaceutical industries. These flasks provide an ideal environment for both adherent and suspension cells, featuring a vented cap that ensures optimal gas exchange, crucial for maintaining healthy cell cultures. The geometric neck design allows easy access for media changes and cell passaging, while the transparent material enables continuous monitoring of cell growth and media conditions. They are made from high-quality polystyrene, offering exceptional chemical resistance and durability, ensuring reliable performance throughout extended cell culture periods.



VISIBLE WRITE-ON SPOT

CLEAR GRADUATIONS



**HYDROPHOBIC
PTFE MEMBRANE**

GEOMETRIC NECK



- ◆ Made of premium polystyrene raw material.
- ◆ Certified free from DNase, RNase and Non-pyrogen.
- ◆ Uniform wall thickness for consistent performance.
- ◆ Clear graduations facilitate easy sample identification.



600mL



250mL



50mL



TC Treated



Non-Treated


3 general sizes of flasks for your choice: 25cm²/75cm²/175cm²

Cell Culture Flask




Cat. No.	Growth Area	Volume	Cap Type	Surface Treatment	Sterile	Packaging
PK00260	25 cm ²	50 mL	Seal Cap	TC Treated	YES	12pieces/bag. 25bags/carton
PK00261	75 cm ²	250 mL	Seal Cap	TC Treated	YES	5pieces/bag. 18bags/carton
PK00262	175 cm ²	600 mL	Seal Cap	TC Treated	YES	5pieces/bag. 10bags/carton
PK00263	25 cm ²	50 mL	Vent Cap	TC Treated	YES	12pieces/bag. 25bags/carton
PK00264	75 cm ²	250 mL	Vent Cap	TC Treated	YES	5pieces/bag. 18bags/carton
PK00265	175 cm ²	600 mL	Vent Cap	TC Treated	YES	5pieces/bag. 10bags/carton
PK00433	25 cm ²	50 mL	Vent Cap	Non-Treated	YES	12pieces/bag. 25bags/carton
PK00434	75 cm ²	250 mL	Vent Cap	Non-Treated	YES	5pieces/bag. 18bags/carton
PK00435	175 cm ²	600 mL	Vent Cap	Non-Treated	YES	5pieces/bag. 10bags/carton
PK00436	25 cm ²	50 mL	Seal Cap	Non-Treated	YES	12pieces/bag. 25bags/carton
PK00437	75 cm ²	250 mL	Seal Cap	Non-Treated	YES	5pieces/bag. 18bags/carton
PK00438	175 cm ²	600 mL	Seal Cap	Non-Treated	YES	5pieces/bag. 10bags/carton

03 CELL CULTURE



CELL CULTURE DISH

PEKYBIO Cell Culture Dishes are primarily used for the cultivation and analysis of cells in controlled conditions, making them vital for both research and diagnostics. With a shallow design and ample surface area, they support easy manipulation and accurate observation of cellular growth. PEKYBIO Cell Culture Dish is made from high-quality polystyrene, offering excellent transparency for clear viewing of cells, and its chemical resistance ensures minimal interference with experimental results, providing a stable environment for a variety of cell culture applications.



**HIGH LIGHT
TRANSMITTANCE**

EASY TO GRIP



- ◆ High clarity, 100 % virgin polystyrene.
- ◆ "TC" refers to vacuum plasma tissue culture surface treatment.
- ◆ Rim on the bottom side of the lid mates with the dish brim for easy and safe stacking.
- ◆ Sterilized by irradiation, SAL= 10^{-6} .




35mm



60mm



100mm




150mm

Cell Culture Dish



Cat. No.	Dimension	Growth Area	Color	Surface Treatment	Sterile	Packaging
PK00266	35Φx12(H)	9.7 cm ²	Clear	TC Treated	YES	20pieces/bag. 25bags/carton
PK00267	60Φx15(H)	22.7 cm ²	Clear	TC Treated	YES	20pieces/bag. 25bags/carton
PK00268	100Φx20(H)	59.4 cm ²	Clear	TC Treated	YES	10pieces/bag. 30bags/carton
PK00269	150Φx25(H)	148.5 cm ²	Clear	TC Treated	YES	10pieces/bag. 8bags/carton
PK00439	35Φx12(H)	9.7 cm ²	Clear	Non-Treated	YES	20pieces/bag. 25bags/carton
PK00440	60Φx15(H)	22.7 cm ²	Clear	Non-Treated	YES	20pieces/bag. 25bags/carton
PK00441	100Φx20(H)	59.4 cm ²	Clear	Non-Treated	YES	10pieces/bag. 30bags/carton
PK00624	150Φx25(H)	148.5 cm ²	Clear	Non-Treated	YES	10pieces/bag. 8bags/carton



PETRI DISH

PEKYBIO Petri Dishes are essential tools in microbiology, cell culture, and biological research. Widely used for the cultivation and observation of bacteria, fungi, and other microorganisms, they provide reliable support for various biological experiments. The transparent material allows for easy observation, with excellent chemical resistance and contamination protection. PEKYBIO Petri Dish is made from high-quality polystyrene, ensuring stability and clarity during use, preventing material-related issues that could affect experimental results.

Designed for everyday laboratory use, the Petri dish meets the needs of bacterial culturing, cell seeding, drug screening, and other applications, ensuring that no external interference affects the experimental process. Ideal for microbiological research and other experiments requiring reliable cultivation conditions.

- ◆ High clarity, 100% virgin polystyrene for clear observation.
- ◆ Two compartments petri dish for dual-sample experiments.
- ◆ Contact petri dish features a grid base for accurate colony counting.
- ◆ Rim design ensures secure stacking and easy handling.
- ◆ Available in multiple specifications to meet diverse experimental needs.
- ◆ Sterilized by irradiation, SAL=10⁻⁶.

03 CELL CULTURE



Vented Lid



PK00938



PK00968



PK00622



PK00705



PK00271



PK01244

Petri Dish



Cat. No.	Dimension	Growth Area	Vents	Material	Compartments	Surface Treatment	Sterile	Packaging
PK00622	60Φ×15(H)mm	23.2 cm ²	3	PS	1	Non-Treated	YES	10pieces/bag. 100bags/carton
PK00705	90Φ×15(H)mm	54 cm ²	4	PS	1	Non-Treated	YES	20pieces/bag. 25bags/carton
PK00271	150Φ×15(H)mm	157.6 cm ²	3	PS	1	Non-Treated	YES	10pieces/bag. 20bags/carton
PK00938	100×100×17(H)mm	84.9 cm ²	3	PS	1	Non-Treated	YES	10pieces/bag. 50bags/carton
PK00968	130×130×17.3(H)mm	136.1 cm ²	4	PS	1	Non-Treated	YES	10pieces/bag. 30bags/carton
PK01244	90Φ×15(H)mm	56.7 cm ²	3	PS	2	Non-Treated	YES	20pieces/bag. 25bags/carton



Contact Petri Dish



Cat. No.	Dimension	Growth Area	Vents	Material	Surface Treatment	Sterile	Packaging
PK01245	55Φ×15(H)mm	25 cm ²	4	PS	Non-Treated	YES	20pieces/bag. 25bags/carton



ERLENMEYER FLASK

PEKYBIO PC/ PETG Erlenmeyer Flasks are ideal for suspension of cell culture, media preparation, mixing and storage. PC/ PETG is robust and break-resistant, increasing personal safety. All sizes are available vented or non-vented. They are non-toxic and environmentally friendly. It is packed individually in a sterile sealable bag. Erlenmeyer flask with vent cap comes with 0.22µm hydrophobic PTFE membrane for continuous gas exchange without contamination.

PC Erlenmeyer Flask



UNIFORM WALL THICKNESS

HIGH TRANSPARENCY



PREMIUM RAW MATERIAL

CLEAR GRADUATIONS



- ◆ PC Bottle & HDPE Cap, high transparency and mechanical strength.
- ◆ Seal cap and vent cap are available.
- ◆ Volume options – 125, 250, 500, 1000mL.
- ◆ Non-Pyrogenic, DNase/RNase free.



Baffled Bottom

03 CELL CULTURE

PC Erlenmeyer Flask- 125, 250, 500, 1000 mL



Cat. No.	Volume	Type of Cap	Color of Cap	Sterile	Type of Bottom	Packaging	
PK01140	125 mL	Vent Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK01141	250 mL	Vent Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK01142	500 mL	Vent Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK01143	1000 mL	Vent Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK01144	125 mL	Seal Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK01145	250 mL	Seal Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK01146	500 mL	Seal Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK01147	1000 mL	Seal Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK01173	125 mL	Vent Cap	Orange	YES	Baffled Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK01174	250 mL	Vent Cap	Orange	YES	Baffled Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK01175	500 mL	Vent Cap	Orange	YES	Baffled Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK01176	1000 mL	Vent Cap	Orange	YES	Baffled Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK01177	125 mL	Seal Cap	Orange	YES	Baffled Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK01178	250 mL	Seal Cap	Orange	YES	Baffled Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK01179	500 mL	Seal Cap	Orange	YES	Baffled Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK01180	1000 mL	Seal Cap	Orange	YES	Baffled Bottom	Individual Packaging	1piece/bag. 25bags/carton

PETG Erlenmeyer Flask

- ◆ PETG Bottle & HDPE Cap.
- ◆ Graduated for quick volume assessment.
- ◆ Vent cap features 0.22µm hydrophobic PTFE membrane for continuous gas exchange.
- ◆ Sterilized by irradiation, SAL=10⁻⁶.
- ◆ Packed in individual bags for easy storage and handling.



PETG Erlenmeyer Flask- 125, 250, 500, 1000 mL



Cat. No.	Volume	Type of Cap	Color of Cap	Sterile	Type of Bottom	Packaging	
PK00272	125 mL	Vent Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK00273	250 mL	Vent Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK00274	500 mL	Vent Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK00275	1000 mL	Vent Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK00276	125 mL	Seal Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK00277	250 mL	Seal Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 50bags/carton
PK00278	500 mL	Seal Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 25bags/carton
PK00279	1000 mL	Seal Cap	Orange	YES	Flat Bottom	Individual Packaging	1piece/bag. 25bags/carton



04 SAMPLE STORAGE

Media Bottle- 54-55

Reagent Bottle- 56-60

04 SAMPLE STORAGE



MEDIA BOTTLE

PEKYBIO PET/ PETG Media Bottles are manufactured in a class 100,000 cleanroom environment. They are ideal for storage and sampling of media, biological and other aqueous solutions. These light-weight bottles are shatter-resistant, sterile and come in a space-saving square design. It can be used for reagent storage but is not suitable for use as a heating vessel. Bottle and cap are normally pre-assembled together.



BAG PACKAGING

PAPER TRAY PACKAGING



INDIVIDUAL PACKAGING

HIGH TRANSPARENCY



- ◆ Clear graduations for approximate volume guidance.
- ◆ Square design allows efficient use of storage space.
- ◆ PET/ PETG bottle with natural HDPE closure, leakproof.
- ◆ Sterilized by irradiation, SAL=10⁻⁶.

PET/ PETG Square Media Bottle with HDPE Closure



PET Square Bottle

PET Square Bottle-50, 60, 125, 250, 500, 1000 mL



Cat. No.	Volume	Type of Bottle	Type of Cap	Pre-Assembled	Sterile	Packaging
PK00384	50 mL	Square	Screw Cap	YES	YES	72pieces/bag. 5bags/carton
PK00281	60 mL	Square	Screw Cap	YES	YES	56pieces/bag. 2bags/carton
PK00282	125 mL	Square	Screw Cap	YES	YES	36pieces/bag. 4bags/carton
PK00283	250 mL	Square	Screw Cap	YES	YES	32pieces/bag. 3bags/carton
PK00284	500 mL	Square	Screw Cap	YES	YES	45pieces/bag. 3bags/carton
PK00285	1000 mL	Square	Screw Cap	YES	YES	32pieces/bag. 2bags/carton

PETG Square Bottle

PETG Square Bottle-30, 60, 125, 250, 500, 1000 mL



Cat. No.	Volume	Type of Bottle	Type of Cap	Pre-Assembled	Sterile	Packaging
PK00926	30 mL	Square	Screw Cap	YES	YES	40pieces/bag. 10bags/carton
PK00929	60 mL	Square	Screw Cap	YES	YES	72pieces/bag. 5bags/carton
PK00532	125 mL	Square	Screw Cap	YES	YES	24pieces/tray. 2trays/bag. 2bags/carton
PK00533	250 mL	Square	Screw Cap	YES	YES	24pieces/tray. 2trays/bag. 2bags/carton
PK00534	500 mL	Square	Screw Cap	YES	YES	24pieces/tray. 2trays/bag. 1bag/carton
PK00535	1000 mL	Square	Screw Cap	YES	YES	12pieces/tray. 2trays/bag. 1bag/carton

04 SAMPLE STORAGE



REAGENT BOTTLE

HDPE Reagent Bottle



PEKYBIO HDPE Bottles are natural and amber, pliable and have excellent strength, including polypropylene caps. They meet the requirements of the most demanding applications in laboratory, such as storing, shipping and packing samples. Amber bottles reduce UV light transmission to protect light-sensitive liquids. HDPE bottle provides excellent chemical resistance, linerless caps provide the leakproof protection, and they are designed for easy and comfortable opening and closing.



Caps in various colors are available.

- ◆ High-density polyethylene bottle with polypropylene seal screw cap.
- ◆ Resistant to most chemicals.
- ◆ Leakproof, linerless caps for secure sealing.
- ◆ Safe, shatter-resistant alternative to glass bottles.
- ◆ Sterile option available.



HDPE Reagent Bottle (Narrow Mouth)



Cat. No.	Volume	Color	Mouth	Pre-Assembled	Dimension	Sterile	Packaging
PK00520	4 mL	Natural	Narrow	NO	38.0 mm x 15.4 mm	NO	200pieces/bag. 10bags/carton
PK00348	8 mL	Natural	Narrow	NO	41.4 mm x 24.6 mm	NO	100pieces/bag. 10bags/carton
PK00349	15 mL	Natural	Narrow	NO	55.3 mm x 24.2 mm	NO	100pieces/bag. 10bags/carton
PK00302	30 mL	Natural	Narrow	NO	58.4 mm x 33.7 mm	NO	100pieces/bag. 10bags/carton
PK00303	60 mL	Natural	Narrow	NO	81.0 mm x 38.6 mm	NO	50pieces/bag. 10bags/carton
PK00304	125 mL	Natural	Narrow	NO	94.2 mm x 49.7 mm	NO	25pieces/bag. 10bags/carton
PK00305	250 mL	Natural	Narrow	NO	127.7 mm x 60.1 mm	NO	20pieces/bag. 10bags/carton
PK00306	500 mL	Natural	Narrow	NO	163.3 mm x 71.7 mm	NO	10pieces/bag. 10bags/carton
PK00307	1000 mL	Natural	Narrow	NO	193.5 mm x 89.8 mm	NO	5pieces/bag. 10bags/carton
PK00522	4 mL	Amber	Narrow	NO	38.0 mm x 15.4 mm	NO	200pieces/bag. 10bags/carton
PK00350	8 mL	Amber	Narrow	NO	41.4 mm x 24.6 mm	NO	100pieces/bag. 10bags/carton
PK00351	15 mL	Amber	Narrow	NO	55.3 mm x 24.2 mm	NO	100pieces/bag. 10bags/carton
PK00308	30 mL	Amber	Narrow	NO	58.4 mm x 33.7 mm	NO	100pieces/bag. 10bags/carton
PK00309	60 mL	Amber	Narrow	NO	81.0 mm x 38.6 mm	NO	50pieces/bag. 10bags/carton
PK00310	125 mL	Amber	Narrow	NO	94.2 mm x 49.7 mm	NO	25pieces/bag. 10bags/carton
PK00311	250 mL	Amber	Narrow	NO	127.7 mm x 60.1 mm	NO	20pieces/bag. 10bags/carton
PK00312	500 mL	Amber	Narrow	NO	163.3 mm x 71.7 mm	NO	10pieces/bag. 10bags/carton
PK00313	1000 mL	Amber	Narrow	NO	193.5 mm x 89.8 mm	NO	5pieces/bag. 10bags/carton

04 SAMPLE STORAGE



HDPE Reagent Bottle (Wide Mouth)



Cat. No.	Volume	Color	Mouth	Pre-Assembled	Dimension	Sterile	Packaging
PK00519	4 mL	Natural	Wide	NO	38.0 mm x 15.4 mm	NO	200pieces/bag. 10bags/carton
PK00286	8 mL	Natural	Wide	NO	41.4 mm x 24.6 mm	NO	100pieces/bag. 10bags/carton
PK00287	15 mL	Natural	Wide	NO	55.3 mm x 24.2 mm	NO	100pieces/bag. 10bags/carton
PK00288	30 mL	Natural	Wide	NO	58.4 mm x 33.7 mm	NO	100pieces/bag. 10bags/carton
PK00289	60 mL	Natural	Wide	NO	81.0 mm x 38.6 mm	NO	50pieces/bag. 10bags/carton
PK01151	100 mL	Natural	Wide	NO	93.5 mm x 45.8 mm	NO	25pieces/bag. 10bags/carton
PK00290	125 mL	Natural	Wide	NO	94.2 mm x 49.7 mm	NO	25pieces/bag. 10bags/carton
PK00291	250 mL	Natural	Wide	NO	127.7 mm x 60.1 mm	NO	20pieces/bag. 10bags/carton
PK00292	500 mL	Natural	Wide	NO	163.3 mm x 71.7 mm	NO	10pieces/bag. 10bags/carton
PK00293	1000 mL	Natural	Wide	NO	193.5 mm x 89.8 mm	NO	5pieces/bag. 10bags/carton
PK00521	4 mL	Amber	Wide	NO	38.0 mm x 15.4 mm	NO	200pieces/bag. 10bags/carton
PK00294	8 mL	Amber	Wide	NO	41.4 mm x 24.6 mm	NO	100pieces/bag. 10bags/carton
PK00295	15 mL	Amber	Wide	NO	55.3 mm x 24.2 mm	NO	100pieces/bag. 10bags/carton
PK00296	30 mL	Amber	Wide	NO	58.4 mm x 33.7 mm	NO	100pieces/bag. 10bags/carton
PK00297	60 mL	Amber	Wide	NO	81.0 mm x 38.6 mm	NO	50pieces/bag. 10bags/carton
PK01152	100 mL	Amber	Wide	NO	93.5 mm x 45.8 mm	NO	25pieces/bag. 10bags/carton
PK00298	125 mL	Amber	Wide	NO	94.2 mm x 49.7 mm	NO	25pieces/bag. 10bags/carton
PK00299	250 mL	Amber	Wide	NO	127.7 mm x 60.1 mm	NO	20pieces/bag. 10bags/carton
PK00300	500 mL	Amber	Wide	NO	163.3 mm x 71.7 mm	NO	10pieces/bag. 10bags/carton
PK00301	1000 mL	Amber	Wide	NO	193.5 mm x 89.8 mm	NO	5pieces/bag. 10bags/carton

PP Reagent Bottle



PEKYBIO PP Reagent Bottles are expertly crafted for laboratory use and offers unparalleled compatibility with a range of lab reagents, including acids, bases, alcohols, aldehydes, esters and aliphatic hydrocarbons.

Our translucent bottles provide exceptional content clarity and can withstand autoclaving when filled or empty. The polypropylene (PP) caps ensure maximum leak protection.

Our bottles are composed of premium polypropylene (PP) material, guaranteeing superior product performance and reliability.



Caps in various colors are available.

- ◆ Wear-resistant, not easily deformed, and highly resistant to pressure.
- ◆ Perfect chemical stability, resistant to 121 °C autoclaving and corrosion.
- ◆ Strong seal to prevent leakage of bottle contents.
- ◆ Multiple specifications: wide /narrow mouth options, translucent /amber colors, and capacities ranging from 4 mL to 1000 mL.
- ◆ Sterile option available.




PP Reagent Bottle (Narrow Mouth)



Cat. No.	Volume	Color	Mouth	Pre-Assembled	Dimension	Sterile	Packaging
PK00447	4 mL	Translucent	Narrow	NO	38.0 mm x 15.4 mm	NO	200pieces/bag. 10bags/carton
PK00448	8 mL	Translucent	Narrow	NO	41.4 mm x 24.6 mm	NO	100pieces/bag. 10bags/carton
PK00449	15 mL	Translucent	Narrow	NO	55.3 mm x 24.2 mm	NO	100pieces/bag. 10bags/carton
PK00450	30 mL	Translucent	Narrow	NO	58.4 mm x 33.7 mm	NO	100pieces/bag. 10bags/carton
PK00451	60 mL	Translucent	Narrow	NO	81.0 mm x 38.6 mm	NO	50pieces/bag. 10bags/carton
PK00452	125 mL	Translucent	Narrow	NO	94.2 mm x 49.7 mm	NO	25pieces/bag. 10bags/carton
PK00453	250 mL	Translucent	Narrow	NO	127.7 mm x 60.1 mm	NO	20pieces/bag. 10bags/carton
PK00454	500 mL	Translucent	Narrow	NO	163.3 mm x 71.7 mm	NO	10pieces/bag. 10bags/carton
PK00455	1000 mL	Translucent	Narrow	NO	193.5 mm x 89.8 mm	NO	5pieces/bag. 10bags/carton
PK00465	4 mL	Amber	Narrow	NO	38.0 mm x 15.4 mm	NO	200pieces/bag. 10bags/carton
PK00466	8 mL	Amber	Narrow	NO	41.4 mm x 24.6 mm	NO	100pieces/bag. 10bags/carton
PK00467	15 mL	Amber	Narrow	NO	55.3 mm x 24.2 mm	NO	100pieces/bag. 10bags/carton
PK00468	30 mL	Amber	Narrow	NO	58.4 mm x 33.7 mm	NO	100pieces/bag. 10bags/carton
PK00469	60 mL	Amber	Narrow	NO	81.0 mm x 38.6 mm	NO	50pieces/bag. 10bags/carton
PK00470	125 mL	Amber	Narrow	NO	94.2 mm x 49.7 mm	NO	25pieces/bag. 10bags/carton
PK00471	250 mL	Amber	Narrow	NO	127.7 mm x 60.1 mm	NO	20pieces/bag. 10bags/carton
PK00472	500 mL	Amber	Narrow	NO	163.3 mm x 71.7 mm	NO	10pieces/bag. 10bags/carton
PK00473	1000 mL	Amber	Narrow	NO	193.5 mm x 89.8 mm	NO	5pieces/bag. 10bags/carton

04 SAMPLE STORAGE



PP Reagent Bottle (Wide Mouth)




Cat. No.	Volume	Color	Mouth	Pre-Assembled	Dimension	Sterile	Packaging
PK00456	4 mL	Translucent	Wide	NO	38.0 mm x 15.4 mm	NO	200pieces/bag. 10bags/carton
PK00457	8 mL	Translucent	Wide	NO	41.4 mm x 24.6 mm	NO	100pieces/bag. 10bags/carton
PK00458	15 mL	Translucent	Wide	NO	55.3 mm x 24.2 mm	NO	100pieces/bag. 10bags/carton
PK00459	30 mL	Translucent	Wide	NO	58.4 mm x 33.7 mm	NO	100pieces/bag. 10bags/carton
PK00460	60 mL	Translucent	Wide	NO	81.0 mm x 38.6 mm	NO	50pieces/bag. 10bags/carton
PK01155	100 mL	Translucent	Wide	NO	93.5 mm x 45.8 mm	NO	25pieces/bag. 10bags/carton
PK00461	125 mL	Translucent	Wide	NO	94.2 mm x 49.7 mm	NO	25pieces/bag. 10bags/carton
PK00462	250 mL	Translucent	Wide	NO	127.7 mm x 60.1 mm	NO	20pieces/bag. 10bags/carton
PK00463	500 mL	Translucent	Wide	NO	163.3 mm x 71.7 mm	NO	10pieces/bag. 10bags/carton
PK00464	1000 mL	Translucent	Wide	NO	193.5 mm x 89.8 mm	NO	5pieces/bag. 10bags/carton
PK00474	4 mL	Amber	Wide	NO	38.0 mm x 15.4 mm	NO	200pieces/bag. 10bags/carton
PK00475	8 mL	Amber	Wide	NO	41.4 mm x 24.6 mm	NO	100pieces/bag. 10bags/carton
PK00476	15 mL	Amber	Wide	NO	55.3 mm x 24.2 mm	NO	100pieces/bag. 10bags/carton
PK00477	30 mL	Amber	Wide	NO	58.4 mm x 33.7 mm	NO	100pieces/bag. 10bags/carton
PK00478	60 mL	Amber	Wide	NO	81.0 mm x 38.6 mm	NO	50pieces/bag. 10bags/carton
PK01156	100 mL	Amber	Wide	NO	93.5 mm x 45.8 mm	NO	25pieces/bag. 10bags/carton
PK00479	125 mL	Amber	Wide	NO	94.2 mm x 49.7 mm	NO	25pieces/bag. 10bags/carton
PK00480	250 mL	Amber	Wide	NO	127.7 mm x 60.1 mm	NO	20pieces/bag. 10bags/carton
PK00481	500 mL	Amber	Wide	NO	163.3 mm x 71.7 mm	NO	10pieces/bag. 10bags/carton
PK00482	1000 mL	Amber	Wide	NO	193.5 mm x 89.8 mm	NO	5pieces/bag. 10bags/carton


Worldwide Distributors Wanted

Ready to Lead Your Market?
Contact Us to Get Partnership Details.





OUR PEKYBIO YOUR LAB SOLUTIONS!



+86 571 8756 6925



+86 571 8756 6925



www.pekybio.com



marketing@pekybio.com



Address.: Building 4, No. 1 Yunhai Road, Lihai Subdistrict
Yuecheng District, Shaoxing City
Zhejiang Province, 312000, China

PEKYBIO logo is the registered trademark of Zhejiang Peki Biotech Co., Ltd.
All rights reserved, including graphics and images.
Copyright 2026 by Zhejiang Peki Biotech Co., Ltd.
The 5th version printed in China 01/2026.

